

Double Pump RF

Medical Vision is introducing the next generation of Double Pump, Double Pump RF, now compatible with tools the surgeon needs during an arthroscopic procedure, Shaver and RF units. This means that suction is simultaneously directed through either the Shaver or RF handpiece.

The compatibility with an RF system is another unique differentiator between Medical Vision® product offering and other available pumps on the market. The other features are:

- Seamlessly integrated in any arthroscopic tower. Now compatible with Shaver systems as well as RF systems.
- Accurate pressure measurement, displaying the true joint pressure through patented algorithms.
- Automated pressure and flow control through the revolutionary SmartVision® technology.
- Easy-to-use Cassette technology safeguards correct mounting of the tubings.
- The Outflow Tracking™ function.



For orthopedic surgeons who expect the highest level of performance and reliability from his Fluid Management system when performing complex arthroscopic interventions, Double Pump RF is the first arthroscopic pump with connectivity to both Shaver and RF systems. Double Pump RF is designed to optimize pressure and flow even in long and challenging procedures.

Double Pump RF

ORDERING INFORMATION:

Double Pump RF	
Item number	Description
1028511	Double Pump RF

Double Pump RF Cassettes		
Item number	Description	Packaging
1002053	Day Cassette*	20 pcs/Box
1003059	RF Patient Cassette	20 pcs/Box
1002052	Single*	100 pcs/Box

*Same as for Double Pump